

Benzophenone-d5 100 µg/mL in Methanol

Art. ID LM24-N-13400-4523D-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzophenone-2,3,4,5,6- d5	[2694-78-2]	100	µg/mL		