

2-Pentanone-1,1,1,3,3-d5 1000 µg/mL in Methanol

Art. ID LM24-N-13600-0929-1000ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Pentanone-1,1,1,3,3-D 5	[24313-49-3]	1000	µg/mL		