

### **Propionic acid-13C3 (C3:0) 100 µg/mL in Acetone**

Art. ID LM24-N-20000-0011C-100AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Propionic acid-13C3	[153474-35-2]	100	µg/mL		