

## **1,1,2,2-Tetrachloroethane-d2 100 µg/mL in Methanol**

Art. ID LM24-N-10900-0711D-100ME2  
Unit 2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1,2,2-Tetrachloroethane-d2	[33685-54-0]	100	µg/mL		