

### **1,2,4-Triazole-d3 100 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-4939D-100AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,4-Triazole-D3	[43088-92-2]	100	µg/mL		