

## **Cyprodinil-d5 10 µg/mL in Acetone**

Art. ID LM24-N-15900-1123D-10AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cyprodinil-(phenyl-d5)	[1773496-67-5]	10	µg/mL		