

**Di-n-propyl phthalate-d4 100 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID LM24-N-13800-0241D-100IO5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-propyl ester D4 (Bis-propyl phthalate)	[358731-29-0]	100	µg/mL		