

### **Irbesartan-d6 10 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-2493D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Irbesartan-d6	[2375621-21-7]	10	µg/mL		