

Isoproturon-d3 10 µg/mL in Acetone

Art. ID LM24-N-15900-5707D-10AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|-------|--------|--------|
| | Isoproturon-d3 | [352438-80-3] | 10 | µg/mL | | |