

## **Ketoprofen-d3 10 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-2494D-10AN10  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ketoprofen D3	[159490-55-8]	10	µg/mL		