

## **Ketoprofen-d3 10 µg/mL in Acetone**

Art. ID LM24-N-18000-2494D-10AC10  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number    | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|---------------|---------------------|-------|--------|--------|
|                  | Ketoprofen D3     | [159490-55-8] | 10                  | µg/mL |        |        |