

Picolinafen-d4 10 µg/mL in Acetonitrile

Art. ID LM24-N-15900-3956D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------------------------|---------------------|-------|--------|--------|
| | Picolinafen-d4 | [137641-05-5]/u nlabelled | 10 | µg/mL | | |