

## **Prometryn-d14 10 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-4091D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Prometryn-D14 (Di-Iso-P ropl-D14)	[1202864-57-0]	10	µg/mL		