

Propanil-d5 100 µg/mL in Acetonitrile

Art. ID LM24-N-15900-4121-100AN1
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Propanil-d5 | [1398065-82-1] | 100 | µg/mL | | |