

Dimethyl phthalate-d4 100 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID LM24-N-13800-0161D-100IO5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-methyl ester D4 (Bis-methyl phthalate D4)	[93951-89-4]	100	µg/mL		