

## **Linuron-d6 10 µg/mL in Acetone**

Art. ID LM24-N-15900-3011D-10A5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Linuron-(methyl-d3, met hoxy-d3)	[1219804-76-8]	10	µg/mL		