

## **Phosmet-oxon-d6 10 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-3923D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosmet-oxon-d6	[n/a]	10	µg/mL		