

Acetylsalicylic acid-d3 10 µg/mL in Acetonitrile

Art. ID LM24-N-18000-2539D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetylsalicylic acid-(methyl-d3)	[921943-73-9]	10	µg/mL		