

## **Salicylic acid-d4 10 µg/mL in Methanol**

Art. ID LM24-N-18000-1041D-10ME10  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Salicylic Acid-d4	[78646-17-0]	10	µg/mL		