

### **Custom Solution - 3 components each at 1 mg/L in 2% HNO<sub>3</sub>**

Art. ID LM24-CUS-70928-100ML  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation: 12 months (originally closed bottle), 6 month after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	1	mg/L		
	Lead (Pb)	[7439-92-1]	1	mg/L		
	Tin (Sn)	[7440-31-5]	1	mg/L		