

Custom Solution - 3 components each at 10 µg/mL in 2% HNO₃

Art. ID LM24-CUS-70930-100ML
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation: 12 months (originally closed bottle), 6 month after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	10	µg/mL		
	Lead (Pb)	[7439-92-1]	10	µg/mL		
	Tin (Sn)	[7440-31-5]	10	µg/mL		