

**Custom Solution - 3 components 1000 mg/L each of in HNO3 2%**

Art. ID LM24-CUS-72031  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

**Description**

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		