

### **Custom Solution - 4 components 5000 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-72035  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### **Description**

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nickel (Ni)	[7440-02-0]	5000	mg/L		
	Copper (Cu)	[7440-50-8]	5000	mg/L		
	Zinc (Zn)	[7440-66-6]	5000	mg/L		
	Lead (Pb)	[7439-92-1]	5000	mg/L		