

Custom Solution - 4 components 5000 mg/L each of in HNO3 5%

Art. ID LM24-CUS-72036
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Manganese (Mn)	[7439-96-5]	5000	mg/L		
	Selenium (Se)	[7782-49-2]	5000	mg/L		
	Chromium (Cr)	[7440-47-3]	5000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	5000	mg/L		