

## **Seronorm™ Immunoassay Liq Low**

Art. ID	SERO-207305
Unit	12 x 3 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Seronorm™ Immunoassay Liq is the "All-in-one" immunoassay control, combining hormones, cardiac and tumour markers. Based on human serum with no preservatives or stabilizers added /// Please note: Analyte composition and analyte concentration may vary with each LOT or instrument. Refer to the package insert of currently available lots for specific analyte and stability claims.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	17-OH-Progesterone					
	25-OH-Vitamin D, total					
	25-OH-Vitamin D3					
	AFP					
	Aldosterone	[52-39-1]				
	AMH					
	Androstenedione	[63-05-8]				
	Anti-TG					
	Anti-TPO					
	beta2-microglobulin					
	beta-hCG, total					
	CA 125					
	CA 15-3					
	CA 19-9					
	CEA					
	CK-MB					
	Cortisol	[50-23-7]				
	C-Peptide					
	DHEA-Sulfate					
	Digoxin	[20830-75-5]				
	EPO					
	Estradiol	[50-28-2]				
	Ferritin					
	Folate					
	FSH					
	hCG, total					
	HE4					

hGH	
Homocysteine	
IgE, total	
Igf-1	[67763-96-6]
IGFBP-3	
Insulin	
LH	
Methylmalonic Acid	
Myoglobin	
NT-proBNP	
Ostase	
Progesterone	[57-83-0]
Prolactin	
Prostatic acid phosphatase	
PSA, complexed	
PSA, free	
PSA, total	
PTH, intact	
SHBG	
T3, free	
T3, total	
T4, free	
T4, total	
TBG	
Testosterone	[58-22-0]
Theophylline	[58-55-9]
THYROGLOBULIN	[9010-34-8]
Troponin I	
Troponin T	
TSH	
Vitamin B12	