

## **Seronorm™ Cardiac Acute Liquid L-2**

Art. ID	SERO-211205
Unit	6 x 2 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Seronorm Cardiac Acute Liq is the most complete acute control on the market, combining the complete panel of cardiac markers and acute markers to rule out life threatening conditions such as sepsis, brain trauma, thrombosis and ectopic pregnancy /// The product is a liquid material based on human serum, no need for reconstitution. Eliminates potential handling errors and is ready to use right after thawing /// Product documentation is provided with assigned values for all analytes on the most frequently used clinical laboratory and point-of-care test systems on the market /// As a result of lack of standardization, the various methods for the markers Troponin I, NT-proBNP and BNP give widely different results. To achieve better control of these critical markers in your laboratory, the low level is available in two variants (L-1A and L-1B) to ensure relevant levels on the various platforms /// Please note: Analyte composition and analyte concentration may vary with each LOT or instrument. Refer to the package insert of currently available lots for specific analyte and stability claims

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	BNP					
	beta-hCG, total					
	CK					
	CK-MB					
	CRP					
	D-dimer					
	Digoxin	[20830-75-5]				
	hCG, total					
	Homocysteine					
	IL-6					
	Myoglobin					
	NT-proBNP					
	Procalcitonin					
	S-100 beta					
	Troponin I					
	Troponin T					