

## **Seronorm™ Immunoprotein Liq L-1**

Art. ID	SERO-210405
Unit	6 x 1 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Protein control, containing more than 20 proteins, e.g. Cystatin C, anti-CCP, RF, CRP and ASL/ASO /// Contains more than 95% human serum with no preservatives or stabilizers added /// Protein documentation with assigned values on nephelometric and turbidimetric methods /// Two clinically relevant levels – each level available separately /// No need for reconstitution, eliminates potential handling errors and ready to use immediately after thawing /// Please note: Analyte composition and analyte concentration may vary with each LOT or instrument. Refer to the package insert of currently available lots for specific analyte and stability claim

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha1-acid-glycoprotei n					
	alpha1-antitrypsin					
	alpha2-macroglobulin					
	Albumin	[9048-46-8]				
	Anti-CCP					
	Apolipoprotein A1					
	ASO					
	beta2-microglobulin					
	C3c					
	C4					
	Ceruloplasmin					
	CRP					
	Cystatin C					
	Digitoxin	[71-63-6]				
	Ferritin					
	Haptoglobin					
	IgA					
	IgE, total					
	IgG					
	IgM					
	Myoglobin					
	Prealbumin					
	Protein, total					
	RF					
	Transferrin					

