

Autonorm Clin Chem Liq L-3

Art. ID	SERO-202215
Unit	12 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Complete chemistry control, including clinical chemistry analytes, lipids, enzymes, proteins, hormones, cardiac markers and therapeutic drugs /// Also includes bicarbonate /// Contains more than 95 % human serum with no preservatives or stabilizers added /// Product documentation with assayed value list for more than 60 analytes including genuine lipids, proteins (e.g. CRP, Apo A1 and Apo B) and therapeutic drugs /// Available in 3 clinically relevant levels. Each level available separately ///Liquid material - no need for reconstitution, eliminates potential handling errors and is ready for use immediately after thawing /// Please note: Analyte composition and analyte concentration may vary with each LOT or instrument. Refer to the package insert of currently available lots for specific analyte and stability claim

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha1-acid-glycoprotei					
	n					
	alpha1-antitrypsin					
	alpha2-macroglobulin					
	ALAT					
	Albumin	[9048-46-8]				
	ALP					
	Amikacin	[37517-28-5]				
	Amylase, pancreas					
	Amylase, total					
	Apolipoprotein A1					
	Apolipoprotein B					
	ASAT					
	ASO					
	Bicarbonat					
	Bile Acid					
	Bilirubin, direct					
	Bilirubin, total					
	C3c					
	C4					
	Calcium (Ca)	[7440-70-2]				
	Carbamazepine	[298-46-4]				
	Ceruloplasmin					
	Chloride (Cl)					
	Cholesterol, HDL					

Cholesterol, LDL	
Cholesterol, total	
Cholinesterase	
CK	
CK-MB	
Copper (Cu)	[7440-50-8]
Creatinine	[60-27-5]
CRP	
Digoxin	[20830-75-5]
Ethanol	[64-17-5]
Ferritin	
Gentamycine	
GGT	
Glucose	
Haptoglobin	
IgA	
IgG	
IgM	
Iron (Fe)	[7439-89-6]
Kappa free light chains	
Kappa/Lambda ratio	
Lactate	
Lambda free light chain	
s	
LDH	
Lipase	
Lithium (Li)	[7439-93-2]
Magnesium (Mg)	[7439-95-4]
Methotrexate	[59-05-2]
Osmolality	
Paracetamol	[103-90-2]
Phenobarbitone	
Phenytoin	[57-41-0]
Phosphorus (P)	[7723-14-0]
Potassium (K)	[7440-09-7]
Prealbumin	
Protein, total	
Quinidine	[56-54-2]

Salicylate	
Sodium (Na)	[7440-23-5]
Theophylline	[58-55-9]
TIBC	
Tobramycin	[32986-56-4]
Transferrin	
Triglyceride	
UIBC	
Urea	[57-13-6]
Urea nitrogen	[57-13-6]
Uric Acid	[69-93-2]
Valproic acid	[99-66-1]
Vancomycin	[1404-90-6]
Zinc (Zn)	[7440-66-6]