

Autonorm™

Art. ID	SERO-103505
Unit	10 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Quality Control covering the most important analytes in the clinical chemistry laboratory /// Animal serum with no preservatives or stabilizers added /// Product documentation with assayed value list for about 30 analytes /// Analyte levels mainly in the normal range /// Pharmaca free /// Freeze-dried material, offering long shelf life /// Please note: Analyte composition and analyte concentration may vary with each LOT or instrument. Refer to the package insert of currently available lots

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ALAT					
	Albumin	[9048-46-8]				
	ALP					
	Amylase, total					
	ASAT					
	Bilirubin, direct					
	Bilirubin, total					
	Calcium (Ca)	[7440-70-2]				
	Chloride (Cl)					
	Cholesterol, HDL					
	Cholesterol, total					
	Cholinesterase					
	CK					
	CK-MB					
	Coppe					
	Creatinine	[60-27-5]				
	GGT					
	GLDH					
	Glucose					
	Iron (Fe)	[7439-89-6]				
	Lactate					
	LDH					
	Lipase					
	Lithium (Li)	[7439-93-2]				
	Magnesium (Mg)	[7439-95-4]				
	Phosphorus (P)	[7723-14-0]				
	Potassium (K)	[7440-09-7]				

Protein, total	
Sodium (Na)	[7440-23-5]
Triglycerides	
UIBC	
Urea	[57-13-6]
Uric Acid	[69-93-2]
Zinc (Zn)	[7440-66-6]