

## **Seronorm™ CRP Liquid L-2**

Art. ID	SERO-213105
Unit	12 x 1 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

CRP Control of superior quality with > 95 % human serum /// Storage temperature at 2-8 °C /// Straight from the fridge to the instrument - no thawing or reconstitution needed /// 1-month stability in opened vials (2-8 °C) /// Three clinically relevant levels with method specific assigned values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	CRP					