

## **Seronorm™ HbA1c Liquid L-1 (40 mmol/mol Hb)**

Art. ID	SERO-209005
Unit	3 x 1 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Human-based quality control with no stabilizers added /// Produced by in-vitro glycation of human blood, this QC material represents an ethical alternative to other materials on the market /// Product documentation with assayed values for immunological and Point of Care (PoC) methods /// Two clinically relevant levels – each available separately /// Liquid material - no need for reconstitution. Eliminates potential handling errors and is ready to use right after thawing /// Good stability in opened vials

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	HbA1c					