

## **Seronorm™ Human**

Art. ID	SERO-200805
Unit	10 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Multi-component assayed control which covers most analytes analysed on clinical chemistry platforms ///Human serum with no preservatives or stabilizers added /// Product documentation with assigned values for more than 45 analytes, including genuine lipids, routine clinical chemistry analytes, enzymes, proteins, cardiac markers and therapeutic drugs ///Two clinically relevant levels – each level available separately ///Freeze-dried material, offering a long shelf life of up to 4 years /// Please refer to the package insert of currently available lots for specific analyte and stability claims

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ALAT					
	Albumin	[9048-46-8]				
	ALP					
	Amylase, pancreas					
	Amylase, total					
	Apolipoprotein A1					
	Apolipoprotein B					
	ASAT					
	Bile					
	Acid					
	Bilirubin, direct					
	Bilirubin, total					
	Calcium (Ca)	[7440-70-2]				
	Chloride (Cl)					
	Cholesterol, HDL					
	Cholesterol, LDL					
	Cholesterol, total					
	Cholinesterase					
	CK					
	CK-MB					
	Copper (Cu)	[7440-50-8]				
	Creatinine	[60-27-5]				
	CRP					
	Digoxin	[20830-75-5]				
	GGT					
	GLDH					

Glucose	
HBDH	
IgA	
IgG	
IgM	
Iron (Fe)	[7439-89-6]
Lactate	
LDH	
Lipase	
Lithium (Li)	[7439-93-2]
Magnesium (Mg)	[7439-95-4]
NEFA	
Osmolality	
Phosphorus (P)	[7723-14-0]
Potassium (K)	[7440-09-7]
Protein, total	
Sodium (Na)	[7440-23-5]
Theophylline	[58-55-9]
Transferrin	
Triglycerides	
UIBC	
Urea	[57-13-6]
Uric	
Acid	
Zinc (Zn)	[7440-66-6]