

Seronorm™ Paediatric

Art. ID	SERO-208005
Unit	6 x 3 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Accuracy control serum specially made for the newborn /// Contains more than 95 % human serum with no preservatives or stabilizers added /// Product documentation including analytical data with documented traceability on all analytes /// Freeze-dried material, offering long shelf life of up to 4 years /// Please refer to the package insert of currently available lots for specific analyte and stability claim

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ALP					
	Bilirubin, direct					
	Bilirubin, indirect					
	Bilirubin, total					
	Calcium (Ca)	[7440-70-2]				
	Chloride (Cl)					
	Glucose					
	Magnesium (Mg)	[7439-95-4]				
	Phosphorus (P)	[7723-14-0]				
	Potassium (K)	[7440-09-7]				
	Sodium (Na)	[7440-23-5]				