

Custom Solution - 8 components 10 mg/l each of in HNO3 2%

Art. ID LM24-CUS-71437-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Selenium (Se)	[7782-49-2]	10	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Thallium (Tl)	[7440-28-0]	10	mg/L		
	Chromium (Cr)	[7440-47-3]	10	mg/L		
	Cadmium (Cd)	[7440-43-9]	10	mg/L		
	Vanadium (V)	[7440-62-2]	10	mg/L		
	Uranium (U)	[7440-61-1]	10	mg/L		