

ARCONIC - Aluminium base, Aluminum, SUS, wrought, typical analysis, Disc 64 mm x 37 mm

Art. ID ALCOA-SQ-19U
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	~0,2	%		
	Iron (Fe)	[7439-89-6]	~0,1	%		
	Copper (Cu)	[7440-50-8]	~0,5	%		
	Magnesium (Mg)	[7439-95-4]	~1	%		
	Chromium (Cr)	[7440-47-3]	~0,1	%		
	Zinc (Zn)	[7440-66-6]	~0,3	%		
	Arsenic (As)	[7440-38-2]	~0,01	%		
	Gallium (Ga)	[7440-55-3]	~0,02	%		
	Phosphorus (P)	[7723-14-0]	~0,007	%		
	Antimony (Sb)	[7440-36-0]	~0,01	%		
	Scandium (Sc)	[7440-20-2]	~0,2	%		
	Vanadium (V)	[7440-62-2]	~0,004	%		
	Zirconium (Zr)	[7440-67-7]	~0,003	%		