

**Tetrafluoro-2-(heptafluoropropoxy) propanoic acid sodium salt GenX (unlabeled), 100 µg/mL in Methanol**

Art. ID CIL-ULM-11542-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                                             | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------------------------------------------|------------|---------------------|-------|--------|--------|
|                  | Tetrafluoro-2-(heptafluoropropoxy) propanoic acid sodium salt | [n/a]      | 100                 | µg/mL |        |        |