

**see new Batch: BAS-SCRM-670/25      Cast Iron, CRM, Ductile (Nodular) Iron, Block 48**  
**mm x 42 mm x 12 mm**

Art. ID                      BAS-SCRM-670/24  
Unit                        block  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,306 ± 0,021	%		
	Cerium (Ce)	[7440-45-1]	0,0128 ± 0,0013	%		
	Chromium (Cr)	[7440-47-3]	0,494 ± 0,015	%		
	Copper (Cu)	[7440-50-8]	1,011 ± 0,009	%		
	Magnesium (Mg)	[7439-95-4]	0,0617 ± 0,0028	%		
	Manganese (Mn)	[7439-96-5]	0,344 ± 0,009	%		
	Molybdenum (Mo)	[7439-98-7]	0,0193 ± 0,0013	%		
	Nickel (Ni)	[7440-02-0]	1,026 ± 0,009	%		
	Silicon (Si)	[7440-21-3]	2,133 ± 0,014	%		
	Titanium (Ti)	[7440-32-6]	0,1021 ± 0,0028	%		
	Vanadium (V)	[7440-62-2]	0,0294 ± 0,0008	%		
	Sulfur (S)	[7704-34-9]	0,0088 ± 0,0010	%		
	Phosphorus (P)	[7723-14-0]	0,0497 ± 0,0019	%		