

## **Certified Reference Material - Color, Water Supply**

Art. ID ERA-661C  
Unit 30 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

### Description

One 30 mL screw-cap bottle yields up to 200 mL after dilution. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit     | Method | Source |
|------------------|-------------------|------------|---------------------|----------|--------|--------|
|                  | Colour            |            | 10,0 - 75,0         | PC units |        |        |