

Certified Reference Material - Color, Water Supply

Art. ID ERA-661C
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

One 30 mL screw-cap bottle yields up to 200 mL after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Colour		10,0 - 75,0	PC units		