

## **Proficiency Testing - Color, Water Pollution**

|                 |   |
|-----------------|---|
| Art. ID         | ERA-882C  |
| Unit            | 30 mL   |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1 |

### Description

One 30 mL screw-cap bottle yields up to 200 mL after dilution. Use with EPA Methods 110.1, 110.2, and 110.3, Standard Methods 2120B, 2120C, 2120E, or other applicable method.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit     | Method | Source |
|------------------|-------------------|------------|---------------------|----------|--------|--------|
|                  | Colour            |            | 10,0 - 75,0         | PC units |        |        |