

## **Certified Reference Material - Color, Water Pollution**

Art. ID	ERA-1070C
Unit	30 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

### Description

One 30 mL screw-cap bottle yields up to 200 mL after dilution. Use with EPA Methods 110.1, 110.2, and 110.3, Standard Methods 2120B, 2120C, 2120E, or other applicable method. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Colour		10,0 - 75,0	PC units		