

QuiKResponse - Color, Water Pollution

Art. ID ERA-1070CQR
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 30 mL screw-cap bottle yields up to 200 mL after dilution. Use with EPA Methods 110.1, 110.2, and 110.3, Standard Methods 2120B, 2120C, 2120E, or other applicable method.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|----------|--------|--------|
| | Colour | | 10,0 - 75,0 | PC units | | |