

## **QuiKResponse - Color, Water Pollution**

Art. ID	ERA-1070CQR
Unit	30 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 30 mL screw-cap bottle yields up to 200 mL after dilution. Use with EPA Methods 110.1, 110.2, and 110.3, Standard Methods 2120B, 2120C, 2120E, or other applicable method.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Colour		10,0 - 75,0	PC units		