

QuiKResponse - Color, Water Supply

Art. ID ERA-661CQR
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 30 mL screw-cap bottle yields up to 200 mL after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Colour		10,0 - 75,0	PC units		