

Custom Solution - 7 components in HNO3 5%

Art. ID LM24-CUS-72084
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 500 | mg/L | | |
| | Boron (B) | [7440-42-8] | 1000 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 1000 | mg/L | | |
| | Phosphorus (P) | [7723-14-0] | 1000 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 1000 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 1000 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 10000 | mg/L | | |