

Custom Solution - 13 components in HNO₃ 1mol/L

Art. ID	LM24-CUS-72089
Unit	1 x 250 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Selenium (Se)	[7782-49-2]	2	mg/L		
	Boron (B)	[7440-42-8]	5	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	5	mg/L		
	Copper (Cu)	[7440-50-8]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Manganese (Mn)	[7439-96-5]	5	mg/L		
	Molybdenum (Mo)	[7439-98-7]	5	mg/L		
	Zinc (Zn)	[7440-66-6]	5	mg/L		
	Sodium (Na)	[7440-23-5]	500	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Sulfur (S)	[7704-34-9]	500	mg/L		
	Cobalt (Co)	[7440-48-4]	5	µg/mL		