

Method 8327 Surrogate Spiking Mixture in Methanol (with 4 molar equiv. NaOH)

Art. ID CIL-ES-5643
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perfluorobutyric acid (PFBA), sodium salt (2,3,4-13C3, 99%)	[2218-54-4]/unlabeled	1000	ng/mL		
	Perfluoropentanoic acid (PFPeA), sodium salt (pentanoyl-13C5, 99%)	[2706-89-0]/unlabeled	1000	ng/mL		
	Perfluorohexanoic acid (PFHxA), sodium salt (13C6, 99%)	[307-24-4]/unlabeled	1000	ng/mL		
	Perfluoroheptanoic acid (PFHpA), sodium salt (heptanoyl-13C7, 99%)	[20109-59-5]/unlabeled	1000	ng/mL		
	Perfluoro-n-octanoic acid (PFOA) (13C8, 99%)	[1350614-84-4]	1000	ng/mL		
	Perfluorononanoic acid (PFNA) (13C9, 99%)	[375-95-1]/unlabeled	1000	ng/mL		
	Perfluorodecanoic acid (PFDA) (13C9, 99%)	[335-76-2]/unlabeled	1000	ng/mL		
	Perfluoroundecanoic acid (PFUA), sodium salt (13C9, 99%)	[60871-96-7]/unlabeled	1000	ng/mL		
	Perfluorododecanoic acid (PFDoA), sodium salt (dodecanoyl-13C12, 99%)	[307-67-5]/unlabeled	1000	ng/mL		
	Perfluorotetradecanoic acid (PFTeDA) (13C2, 99%)	[376-06-7]/unlabeled	1000	ng/mL		
	Perfluorobutanesulfonate (PFBS), potassium salt (13C4, 99%)	[29420-49-3]/unlabeled	1000	ng/mL		
	Potassium perfluoro-1-h	[3871-99-6]/unlabeled	1000	ng/mL		

hexanesulfonate (PFHxS) (13C6, 99%)	labeled			
Potassium perfluoro-1-octanesulfonate (PFOS) (13C8, 99%)	[2795-39-3]/unlabeled	1000		ng/mL
1H,1H,2H,2H-Perfluorohexanesulfonate (4:2 FTS), sodium salt (13C2, 99%, D4, 98%)	[27619-93-8]/unlabeled	1000		ng/mL
1H,1H,2H,2H-Perfluorooctanesulfonate (6:2 FTS), sodium salt (13C2, 99%, D4, 98%)	[27619-94-9]/unlabeled	1000		ng/mL
1H,1H,2H,2H-Perfluorodecane sulfonate (8:2 FTS), sodium salt (13C2, 99%, D4, 98%)	[27619-96-1]/unlabeled	1000		ng/mL
Perfluorooctanesulfonamide (PFOSA) (13C8, 99%)	[754-91-6]/unlabeled	1000		ng/mL
N-Methylperfluorooctanesulfonamidoacetic acid (N-MeFOSAA) (N-methyl-D3, 98%)	[2355-31-9]/unlabeled	1000		ng/mL
N-Ethylperfluorooctanesulfonamidoacetic acid (N-EtFOSAA) (N-ethyl-D5, 98%)	[2991-50-6]/unlabeled	1000		ng/mL