

**Custom Solution - 5 components 0,1 mg/L each of in HNO3 2%**

Art. ID LM24-CUS-72105  
Unit 2 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

**Description**

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	0,1	mg/L		
	Nickel (Ni)	[7440-02-0]	0,1	mg/L		
	Chromium (Cr)	[7440-47-3]	0,1	mg/L		
	Vanadium (V)	[7440-62-2]	0,1	mg/L		
	Cadmium (Cd)	[7440-43-9]	0,1	mg/L		