

Custom Solution - 4 components 100 mg/L each of in Methanol

Art. ID	LM24-CUS-72107
Unit	3 x 1 mL ampoule
Deliverydetails	De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
technical	Phthalic acid, bis-iso-octyl ester (Bis-iso-octyl phthalate)	[27554-26-3]	100	mg/L		
	Phthalic acid, bis-C6-C8-branched alkyl esters C7-rich	[71888-89-6]	100	mg/L		
mixture of branched chains in isomers	Di-iso-nonyl phthalate	[28553-12-0]	100	mg/L		
	Diisodecyl phthalate (mixture of isomers)	[26761-40-0]	100	mg/L		