

Aluminium base purity 99.999%, RM, cast, Disc 60 mm x 20 mm, Category C

Art. ID SPE-R-11
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	Base			
	Antimony (Sb)	[7440-36-0]	0,000015	%		
	Arsenic (As)	[7440-38-2]	<0,00008	%		
	Barium (Ba)	[7440-39-3]	0,000011	%		
	Beryllium (Be)	[7440-41-7]	<0,000005	%		
	Bismuth (Bi)	[7440-69-9]	0,000001	%		
	Boron (B)	[7440-42-8]	<0,00003	%		
	Cadmium (Cd)	[7440-43-9]	0,0000005	%		
	Calcium (Ca)	[7440-70-2]	0,000015	%		
	Cerium (Ce)	[7440-45-1]	0,000026	%		
	Chromium (Cr)	[7440-47-3]	0,0000017	%		
	Cobalt (Co)	[7440-48-4]	0,0000003	%		
	Copper (Cu)	[7440-50-8]	0,00004	%		
	Gallium (Ga)	[7440-55-3]	0,0000046	%		
	Hafnium (Hf)	[7440-58-6]	0,0000029	%		
	Indium (In)	[7440-74-6]	<0,00001	%		
	Iron (Fe)	[7439-89-6]	0,000072	%		
	Lanthanum (La)	[7439-91-0]	0,0000059	%		
	Lead (Pb)	[7439-92-1]	0,0000011	%		
	Lithium (Li)	[7439-93-2]	<0,00002	%		
	Magnesium (Mg)	[7439-95-4]	0,00009	%		
	Manganese (Mn)	[7439-96-5]	0,0000015	%		
	Molybdenum (Mo)	[7439-98-7]	0,000003	%		
	Nickel (Ni)	[7440-02-0]	0,000001	%		
	Phosphorus (P)	[7723-14-0]	0,0002	%		
	Potassium (K)	[7440-09-7]	<0,00009	%		
	Samarium (Sm)	[7440-19-9]	0,0000084	%		
	Scandium (Sc)	[7440-20-2]	0,000002	%		
	Selenium (Se)	[7782-49-2]	<0,000005	%		
	Silicon (Si)	[7440-21-3]	<0,0001	%		
	Silver (Ag)	[7440-22-4]	0,0000001	%		

Sodium (Na)	[7440-23-5]	<0,00002	%
Strontium (Sr)	[7440-24-6]	0,000008	%
Terbium (Tb)	[7440-27-9]	0,0000008	%
Thorium (Th)	[7440-29-1]	0,0000043	%
Tin (Sn)	[7440-31-5]	<0,00002	%
Titanium (Ti)	[7440-32-6]	<0,00004	%
Uranium (U)	[7440-61-1]	0,0000163	%
Vanadium (V)	[7440-62-2]	<0,00001	%
Zinc (Zn)	[7440-66-6]	<0,00002	%
Zirconium (Zr)	[7440-67-7]	<0,00005	%